KPR Koester Performance Research

Better Computer Access Solutions through the Use of Evidence

Heidi Koester, Ph.D.

hhk@kpronline.com Koester Performance Research Ann Arbor, MI

www.kpronline.com

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Researc

Overview

- Brief Introduction
 - Computer Access Intervention Process
 - · Clear Evidence, Better Solutions
- Use of Evidence in Computer Access
 - · Pointing Devices
 - Text Entry
 - Switch Use

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Access Intervention Process

- Goal: to find the best computer access solution for an individual's needs
- · Determine client needs and goals
- · Assess characteristics of client, environment, tasks
- · Compare possible solutions for input & output
- Recommend & implement particular solution
- · Configure for user's needs
- Measure outcomes

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Clear Evidence, Better Solutions

- Make decisions based on clear evidence related to the client
 - · External or field evidence
 - · Knowledge & skills of the providers
 - · Individual evidence
 - · Discussion with client
 - Functional skills assessment
 - Trials with potential interventions

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

NPR Koester Performance Research

Clear Evidence – example uses

- · Setting goals: therapy, IEP, etc.
 - · Identifying needs
 - · Justifying areas of work
- · Funding support
- · Choosing methods and techniques
- Tracking progress
- · Measuring outcomes

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

Types of Evidence

- Observation
- User Feedback
- Measurement



MPR Koester Performance Research

Measurement Focus

- Assess user abilities that are most relevant to computer usage
- Speed and accuracy during performance of specific fundamental tasks
 - Text entry
 - · Target selection
 - · Target dragging
 - · Switch press

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Resear

Computer-based Measurement Tools

- Present tasks that are repeatable, relevant, and realistic
- Efficient data collection and report generation
- Ideally get the information you need, in less time!
- Focus on KPR software, but will mention other options

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Pointing Device Use

- · Select best pointing method
 - Mice, trackballs, trackpads, head-controlled mice, keyboard-based approaches, etc.
- · Configure for user's needs
 - · Location, splinting, device behavior, etc.
- Follow-along
- · Very little external evidence available

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Skills & Measures for Pointing

- Target Acquisition
 - Dwell, click, double-click
- Dragging
- Menu Selection
- Speed & accuracy
- Complements observations & feedback

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Gathering Pointing Device Evidence

- · Available, but not recommended:
 - · Custom Solutions website
 - · Reach assessment tests
- Better, research-based options (IMHO):
 - · Compass software
 - · Pointing Wizard

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research


Compass Software

- Aim, Drag, and Menu tests specifically for pointing device skills
- Assess on-screen keyboard use with text entry tests
- · Setup is highly customizable, if desired
- · Validity has been demonstrated
- · Compatible with alternative inputs and outputs
- Reports and data are stored for easy review and retrieval

Compass Example

- Selecting the best pointing device for an individual with upper extremity impairments
- Young adult with CP performed Compass Aim tests with three different pointing devices
- Test set-up was identical for each device

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN



Compass Example: Results


	Trial Time (sec)	Entries
Mouse	2.6	1.4
Trackpad	4.9	1.3
Trackball	5.4	1.3


- · View PDF Report...
- · Control looked similar, qualitatively
- · But performance was much faster with the mouse
- Provides team (including user) with means of making an informed decision

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

Pointing Wizard Example

- · Adjusting Windows settings to meet user's needs
- · Pointer speed: gain of the pointing device
- Double-click settings:
 - · Double-click time
 - · Double-click distance
- Object sizes: menus, caption buttons, scrollbar, taskbar





KPR Koester Performance Research

Pointing Summary

- · Clear evidence, Better solutions
- Compass: for assessment, more formal evidence-gathering
- Pointing Wizard: tool for end user to leverage Windows settings effectively

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Text Entry

- Using keyboard or alternatives such as on-screen keyboard, speech recognitions, scanning
- · Some external evidence exists
- Gather internal evidence about the user:
 - · Choose the best text entry method
 - · Configure the text entry method
 - · Track progress

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Skills & Measures for Text Entry

- Reliable and efficient access to all characters and functions
- · Words (or characters) per minute
 - · overall and with errors removed
- Errors
- Complements observations and feedback

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

External Evidence for Typing

Input Method (wpm)
Slow Average Fast
UE Typing, touch 20 30 60

Input Method	Text Entry Performance (wpm)		
mpat motriou	Slow	Average	Fast
UE Typing, touch	20	30	60
Speech Recognition	< 10	16	30
UE Typing, C5-6	4	16	25
Mouthstick Typing	< 5	10 – 12	20
Morse Code	< 5	8	12
On-Screen Keyboard	< 5	6 – 8	10
Row-Col Letter Scanning	< 1		8
UE Typing, MS			
Tongue-touch Keypad	< 1		

How would you use this sort of information?

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Gathering Text Entry Evidence


- · List of tools:
 - · "Standard" typing tests
 - TextTest and StreamAnalyzer
 - · Compass software
 - Keyboard Wizard

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN


OKPR Koester Performance Research


Free Typing Tests


- TypingTest.com
- · speedtest.10-fast-fingers.com
- · Typical limitations:
 - · Distracting screen
 - One size fits all
 - · Can't save results
 - Psychometric validity unknown






















Text Entry Example

• 68 y/o woman with multiple sclerosis

• Runs through Keyboard Wizard (demo)

• Adjusting auto-repeat setting:

• Improved typing speed 50% (from 2.2 to 3.2 wpm)

• Reduced errors 32 pp (from 60% errors to 28%)

Text Entry Summary

Clear evidence, Better solutions
Use external evidence for perspective
Free typing tests OK with caveats

Compass: for assessment, more formal evidence-gathering

Keyboard Wizard: tool for end user to leverage Windows settings effectively

Switch Use

• Select best switch
• Location
• Activation method

• Configure scanning software
• Scan rate

• Compass or SSPT for gathering evidence


KPR Koester Performance Research


Skills & Measures for Switch Use

- · Reliable and efficient switch press
- · Press time
- · Release time
- Errors (extra switch hits)
- Complements observations and feedback

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN









- High school student with cerebral palsy
- No reliable use of a single switch to access a computer or other devices
- Use Compass to compare several switch sites





Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MM

Compass Switch Press Results

	Correct Trials	Avg. Trial Time (s)	Avg. Press Time (s)	Avg. Release Time (s)
Head Right	5/5	8.5	7.8	0.64
Head Left	5/5	9.0	8.5	0.57
Head Posterior	4/5	17.4	11.1	6.32
Right Hand	2/5	19.4	15.6	3.8
Right Finger	4/5	11.2	8.6	2.5

How would you use this information?

MPR Koester Performance Research

Choosing an appropriate scan rate

- Simple estimate:
 - Appropriate scan rate = pressTime * 1.5
- Establishes a scan rate that is generally neither too fast nor too low, for an experienced scanner
- Example: pressTime = 2.5 seconds
 - Scan rate = 2.5 * 1.5 = 3.8 sec

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Besearch

Switch Use Summary

- · Clear evidence, Better solutions
- Measure switch press time for:
 - · Finding best switch site
 - · Establishing an appropriate scan rate

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Conclusions

- Gathering clear evidence can enhance decision-making and document success of an intervention
- Clients themselves often really appreciate seeing data about their performance
- Give one of these tools a try, if you haven't already

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

KPR Software Details

- · Resources at KPR website, kpronline.com
 - Compass, Win/Mac, \$179
 - · Wizards, Win, \$16.95 each
 - Compass + Wizards bundle, \$195
 - · Free trial of Compass available
 - Demonstration videos
- Compass also available through Infogrip, EnableMart, and Technology for Education

Closing the Gap Conference | October 12-14, 2011 | Minneapolis, MN

KPR Koester Performance Research

Final Words

- KPR research & development is supported by:
 - · National Institutes of Health
 - Paralyzed Veterans of America Research Foundation
- · Thanks for being here today!
- Contact: Heidi Koester hhk@kpronline.com